

MATTRESS UNPACKING

Unroll, allow to recover and inspect before use.

01

Check the packaging for damage.

Confirm the model and dimensions and ensure that the packaging is undamaged.

02

Carefully open the packaging.

Remove the wrapping carefully to protect the mattress.

03

Never use knives or sharp objects directly on the mattress.

Such objects may damage the cover and compromise the product's integrity.

04

Place the mattress on a clean, flat surface.

Ensure the surface is dry and free from debris.

05

Unroll the mattress completely.

Unroll carefully without excessive pulling, forced bending or dragging across abrasive surfaces.

06

Allow the foam to fully recover its shape.

The mattress may be used before 48 hours, but 48 hours are required for it to fully recover its shape.

07

Visually inspect the mattress before installation.

Check for visible damage and ensure that the foam has recovered evenly.





4. GOOD PRACTICES FOR USE

-  Use only on compatible beds.
-  Keep the mattress fully supported by the bed base.
-  Always position the foam cuts facing down.
-  Do not bend the mattress excessively.
-  Do not drag the mattress over abrasive surfaces.
-  Do not puncture or cut the cover.
-  Regularly inspect the integrity of the cover and zipper.
-  Remove any damaged mattress from use.



5. CLEANING AND DISINFECTION

- 1 Clean with water and mild soap.
- 2 Use only disinfectants compatible with polyurethane coatings and follow the recommended concentration.
- 3 Always rinse the surface with clean water after disinfection.
- 4 Allow to dry completely before reuse.

Recommended	Not recommended
 Water and mild soap	 Abrasive products
 Soft cloth	 Metal brushes
 Compatible disinfectants	 Steam
 Rinse after disinfection	 Iron
 Dry completely	 Incompatible solvents